

3M™ Scotchcast™ Conformable Roll Splint

Commonly Asked Questions

Question:

What is the advantage of using fiberglass splints vs. plaster?

Answer:

Fiberglass is lightweight and has a lower profile; thus, patient is able to lead a more normal lifestyle. Fiberglass is stronger, there is less splint breakdown, and it is less messy for the clinician.

Question:

How does the Scotchcast splint set time compare to what I'm using now?

Answer:

You control the working time based on your application and needs.

Question:

What education and training is available for your splints?

Answer:

3M offers workshops on both splinting and casting.